




SPC: Datapaq Value (Low 12:00.00 at 190.0°C, Mid 10:00.00 at 200.0°C, High 09:00.00 at 210.0°C) [Full Zoom]



●	#2
▲	#3
▼	#4
+	#5
×	#6

Probe	Maximum	Minimum	Mean	Standard Deviation	Cp		Cpk		Pp		Ppk		Within Spec.
					Per Probe	Overall	Per Probe	Overall	Per Probe	Overall	Per Probe	Overall	
#2 (°C)	111	71	88	8.64	2.11	1.40	0.86	1.07	1.74	1.01	0.71	0.77	X
#3 (°C)	134	86	101	11.68	1.78		1.21		1.28		0.87		X
#4 (°C)	143	93	109	12.21	1.80		1.55		1.23		1.06		X
#5 (°C)	153	96	113	12.81	1.88		1.79		1.17		1.11		X
#6 (°C)	153	94	112	13.84	1.74		1.62		1.08		1.01		X

Outer Upper Limit	160
Inner Upper Limit	120
Target Value	100
Inner Lower Limit	90
Outer Lower Limit	70



Files used for SPC

Name	Directory	Created
Oven SPC Run 1.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	01/11/1993 13:28:00
Oven SPC Run 2.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	02/11/1993 07:28:00
Oven SPC Run 3.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	03/11/1993 07:28:00
Oven SPC Run 4.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	04/11/1993 08:29:00
Oven SPC Run 5.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	05/11/1993 07:29:00
Oven SPC Run 6.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	08/11/1993 07:34:00
Oven SPC Run 7.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	09/11/1993 07:28:00
Oven SPC Run 8.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	10/11/1993 07:32:00
Oven SPC Run 9.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	12/11/1993 07:30:00
Oven SPC Run 10.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	15/11/1993 07:28:00
Oven SPC Run 11.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	16/11/1993 07:28:00
Oven SPC Run 12.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	17/11/1993 07:36:00
Oven SPC Run 13.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	18/11/1993 07:28:00
Oven SPC Run 14.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	19/11/1993 07:30:00
Oven SPC Run 15.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	22/11/1993 08:03:00
Oven SPC Run 16.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	24/11/1993 07:11:00
Oven SPC Run 17.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	25/11/1993 07:21:00
Oven SPC Run 18.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	26/11/1993 07:11:00
Oven SPC Run 19.paq	C:\Datapaq\Insight\OvenProf V7\Paqfiles\	29/11/1993 07:32:00